
# Microarray analysis of alternative transcription






Malachi Griffith
Genome Sciences Centre


GSC Science Retreat - Nov. 2 2006

## Background and rationale



- Alternative transcription generates multiple isoforms from most human loci
  - ~75% of human genes are alternatively spliced (Johnson 2003)
- Specific isoforms may represent useful therapeutic targets or diagnostic markers
- Recent developments in microarray technology allow the efficient detection of specific isoforms
  - A) Affymetrix exon microarrays
  - B) Custom splicing microarrays





Therapeutic antibody specific to CD44 alternate exon 6 (Venables, 2006)

2

## A) Description of Affymetrix exon arrays




- Array design (Human exon 1.0 ST)
  - $-\sim 150,000 loci$
  - $-\sim$ 5.5 million probes (25-mer oligos)
  - $-\sim$ 1.4 million exons (4 probes per exon).
    - ~60% of these exons are 'speculative'
- Samples
  - 1-5 μg total RNA input
- Analysis
  - Identify differential transcription events corresponding to one or more exons

© British Columbia Cancer Agency Branch




# Affymetrix exon array strategy





# 5-FU resistance is associated with overexpression of a short HHIP isoform





HHIP – Hedgehog interacting protein, a HH antagonist

# B) Description of ALEXA platform




- Array design
  - ALEXA platform (Perl, R, mySQL)
    - Works on any EnsEMBL genome
  - Extract probes corresponding to all exons, exon-exon junctions, exon boundaries and introns
    - Vary probe length to achieve target Tm (e.g. 36 bp +/- 10 bp)
  - Filter probes for specificity and thermodynamic properties
  - Select ~400,000 probes and submit to NimbleGen
- Samples
  - 4 μg mRNA (polyA+) input
- Analysis
  - In addition to identifying differential exons, the connections and boundaries of exons are interrogated

© British Columbia Cancer Agency Branch




# ALEXA splicing array strategy





## Affymetrix versus custom ALEXA arrays





## Conclusions



- Preliminary experiments with Affymetrix exon arrays seem promising
  - Identify differential expression of novel and known isoforms. 35 events with 4-fold change or greater
    - 25 affect ORF
    - 9 in 5' UTR and
    - 1 in 3' UTR
  - Low cost, low sample requirements, fast turn-over
- Proof-of-principle custom ALEXA experiments with NimbleGen arrays are also encouraging
  - Moderate cost, high sample requirements, moderate turnover, customizable, open source, species generic
- Experiments to directly compare platforms underway

9

# Acknowledgements



Marco Marra Joseph Connors Stephane Flibotte Steve Jones Gregg Morin

## **Prostate Cancer Group**

- Marianne Sadar
- Tammy Romanuik

#### **Bioinformatics**

- Obi Griffith, Anca Petrescu

## **MGC Group**

- Martin Hirst, Thomas Zeng

#### **GI Cancer Group**

- Isabella Tai, Michelle Tang

#### **Affymetrix Group**

- Agnes Baross, Susanna Chan, Jennifer Asano

#### **Proteomics Group**

- Gregg Morin, Grace Cheng





